



Out-of-the-box label printing: **FLUICS PRINT** for **eLabNext**



eLabNext
EPPENDORF GROUP

FLUICS
P R I N T




Jump-start your label printing!


Print your first labels in only a few minutes. No configuration or coding of templates required.



- 1.** Connect printer to Internet (LAN cable)



- 2.** Activate the FLUICS PRINT add-on from the eLabNext Marketplace



- 3.** Print labels from any PC, Mac, or tablet



What's in the box:

High resolution printer with a small footprint

Pixel perfect QR code printing (300 dpi) with a compact 9.7x6.9x5.3" (245x175x136 mm) footprint. No additional computer required in the lab.



Activation key

Connect your eLabNext account through the cloud to the printer. Print labels from any PC, Mac or tablet.

Best-in-class cryo labels

1,000 labels 1.0x0.5" (25.4x12.7mm), waterproof, cryo from -196°C to +110°C, ultra-permanent, blank, yellow (Type 1). Including black print ribbon and cloud print service.



Expert label designs

Two label designs with tested print settings for optimal scan- and readability (incl. QR code, eLab ID, sample name).

For eLabNext Customers only:

Additional custom label templates free of charge in the first year.



Premium support

Online troubleshooting and 2 years warranty on hardware.

One-time investment. No subscription. Only 1,220 EUR or 1,560 USD.

Additional labels starting at 0.19 EUR or 0.24 USD/label

Order now:

go.fluics.com/elab (EU and ROW) or

go.fluics.com/elab-us (distributor for USA and CAN)



© **FLUICS GmbH**
Fürstenstr. 15
80333 München
Germany
print@fluics.com
www.fluics.com